

Eppendorf Advantage

Advance your research through insights,
solutions, optimal outcomes, and
special offers.

JANUARY 1 – JUNE 30, 2025



Save on Consumables, Centrifuges, Thermal Cyclers, Pipette Bundles, and More!

Ordering information

Description	Part Num.
Get Up to 65% Off on Eppendorf Sample Handling Packages	
2 x twin.tec® <i>real-time</i> PCR Plate 96, unskirted, low profile, white, 20 pcs.	E003013270P
2 x twin.tec® <i>real-time</i> PCR Plate 96, semi-skirted, high profile, blue, 25 pcs.	E951022043P
2 x twin.tec® <i>real-time</i> PCR Plate 96, skirted, low profile, white, 25 pcs.	E951022015P
Masterclear® <i>real-time</i> PCR Film (self-adhesive), 100 pieces	14-285-820P
1 x Mastercycler® X50s, 110-230V 96-well thermoblock, silver with control panel EU and US plug and 2 x Mastercycler® X50i, 110-230V 96-well thermoblock, silver without control panel EU and US plug	E631100001P
1 x Mastercycler® X50a, 110-230V 96-well thermoblock, aluminum with control panel EU and US plug and 2 x Mastercycler® X50l, 110-230V 96-well thermoblock, aluminum without control panel EU and US plug	E631300018P
1 x Mastercycler® X50h, 110-230V 384-well thermoblock, aluminum high pressure lid with control panel, EU and US plug and 2 x Mastercycler® X50t, 110-230V 384-well thermoblock, aluminum high pressure lid without control panel, EU and US plug	E631600019P
Eppendorf ThermoMixer® C, 100 – 130 V/50 – 60 Hz; Eppendorf SmartBlock 1.5 mL, incl. Transfer Rack 1.5/2.0 mL	05-412-503P
Get Up to 45% Off on Eppendorf Liquid Handling Packages	
Xplorer® plus, 8-channel, 5 – 100 µL, yellow, 100 – 240 V ±10 %, 50 – 60 Hz; Xplorer®, 8-channel, 15 – 300 µL, orange, 100 – 240 V ±10 %, 50 – 60 Hz; Xplorer® plus, 8-channel, 50 – 1,200 µL, green, 100 – 240 V ±10 %, 50 – 60 Hz	14-286-558PM
Eppendorf Research® plus, 1-channel, 20 – 200 µL, yellow, variable volume; Eppendorf Research® plus, 1-channel, 0.1 – 2.5 µL, dark gray, variable volume; Eppendorf Research® plus, 1-channel, 100 – 1,000 µL, blue, variable volume; Eppendorf Research® plus, 1-channel, 2 – 20 µL, yellow, variable volume incl. epT.I.P.S.® Box 2.0 with 96 pipette tips; Repeater® M4 Starter Kit, 1-channel, 1 µL – 10 mL, incl. Combitips® advanced Box 2.0, Combitips® advanced assortment pack (1 tip of each size) and Pipette Holder 2	14-287-150PM
Get Up to 55% Off Eppendorf Centrifuge and Shaker Packages	
Microcentrifuge 5425, non-refrigerated, knob control, 120V/50 - 60 Hz, Rotor FA-24x2; includes 3x Eppendorf Research® plus (one-channel) pipettes with 3x epTIPS Box 2.0 (x96 pipette tips)	13-864-454PM1
Vacufuge® plus system with integrated diaphragm vacuum pump, incl. integrated pump, 120 V/50 – 60 Hz (US), w/ Rotor F-45-48-11 & Rotor F-35-6-30; Eppendorf Conical Tubes SnapTec® 50	05-414-177PM
Centrifuge 5430 R, rotary knobs, refrigerated, 120 V/50 – 60 Hz (US); aerosol-tight, incl. rotor lid; incl. wind shield and rotor lid, with 2 buckets; Eppendorf Xplorer®, 8-channel, 15 – 300 µL, orange, 100 – 240 V ±10 %, 50 – 60 Hz	13-690-013PM1
Centrifuge 5920 R, keypad, refrigerated, 120 V/50 – 60 Hz (US), includes Rotor S-4xUniversal-Large, 2 adapters, 2 caps	13-864-306PM1
3x Innova® S44i refrigerated, stackable up to 3-fold, 120 V ± 10 %/50/60 Hz (US), orbit 2.5 cm (1 in); 3x Universal platform; Triple stacking kit	05-400-165PM
Centrifuge 5430 R, keypad, refrigerated, 120 V/50 – 60 Hz with Rotor FA-45-30-11 aerosol-tight rotor, incl. rotor lid; Rotor A-2-MTP incl. wind shield and rotor lid, w/ 2 buckets; Eppendorf Xplorer®, 8-channel, 15 – 300 µL, orange, 100 – 240 V ±10 %, 50 – 60 Hz	13-690-008PM1
Innova 42R with 3/4" orbit, 120V/60Hz; Stacking kit; Universal Platform 18" x 18"	05-400-157PM
Innova 42R with 1" orbit, 120V/60Hz; Stacking kit; Universal Platform 18" x 18"	05-400-162PM

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Call customer service: 1-800-234-7437